



Specifications

Product Code	Diameter	Length	Direction of View
Nasal Endoscopes			
GME0040175	4mm	175mm	0°
GME3040175	4mm	175mm	30°
GME4540175	4mm	175mm	45°
GME7040175	4mm	175mm	70°
Nasal Endoscopes (Paediatric)			
GME0027175	2.7mm	175mm	0°
GME3027175	2.7mm	175mm	30°
Laryngoscopes			
GME9080175	8mm	175mm	90°
GME9060175	6mm	175mm	90°
GME7080175	8mm	175mm	70°
GME7060175	6mm	175mm	70°
Oto-Endoscopes / Tympanoscopes			
GME0040100	4mm	100mm	0°
GME3040100	4mm	100mm	30°
GME0027100	2.7mm	100mm	0°
GME3027100	2.7mm	100mm	30°
GME004050	4mm	50mm	0°
GME304050	4mm	50mm	30°
GME002750	2.7mm	50mm	0°
GME302750	2.7mm	50mm	30°

Features

HD Quality

Wide-Angle

German Optics

Sapphire Scratch Resistant Lens

* Please contact us for more information.

* With a policy of continuous development, we reserve the right to change design and specifications without prior notice.

