


BASIC PLUS

- Compact Ergonomic Design
- Feature Loaded



GEM SURGITECH

refined engineering



ergonomic design

Endoview Endoscopy Camera
LED Monitor with Adjustable Stand

In-Built Suction Machine

- Twin Jars
- Pressure Gauge & Release Valve

High Performance LED Penlight
for Endoscope Illumination

Thermal Cautery

Sliding Drawers
2 Large Size



—○ Standard colors



French Blue
Barn Red



Rich Black
Gold

Technical Details

Suction:	Max Vacuum: 710 $\pm$ 10mm Hg regulable through control knob Suction Capacity: 40 ltrs. / Min Vacuum Gauge: Bourden Type 3" 0-760mm Hg caliberation Motor: 1/4 hp motor Jars: Twin Polycarbonate Jars 2+2 ltrs.
Light Source	Illumination: Halogen Power: 150W Color Temp: 3200K
Anti Fog Device:	Airflow Rate: 40 cu ft./min Air Temperature: 55° C
Endoscopy Camera: Endoview	Image Sensor: Sony Image Sensor 1/3" Signal: PAL Cable Length: 3 metre.
Penlight:	Power: 10W Illumination: LED Color Temp: 6000K
Power Consumption:	500-1000W
Power Supply:	AC 220V-240V, 50/60Hz
Dimensions:	600 x 580 x 1335mm (L x W x H)
Weight:	85kg. approx.

* Please contact us for more information.

* With a policy of continuous development, we reserve the right to change design and specifications without prior notice.



Contact

Address: 155, Preet Nagar, Ambala Cantt., Haryana (133001) India

Contact Info: Ph No. +91-999625-6235, +91-98174-19409

Email: info@gemsurgitech.com, gemsurgitech@rediffmail.com

Website: www.gemsurgitech.com

(An ISO 9001:2015, ISO 13485:2016 Certified Company)

Visit our website:

www.gemsurgitech.com