

Standard Features

- Built-in Twin Halogen Light Source 15v/150w with Brightness Control
- Built-in Suction Machine with Pressure Gauge (Twin Jars)
- Slim LED X-Ray View Box with Adjustable Brightness
- Thermal Cautey
- Anti Fog Device/Mirror Warmer
- 10W High Performance LED Penlight for Endoscope Illumination.
- x2 Stainless Steel Instruments Trays with Cover
- x2 Medicine Bottles with Stand.
- x2 Sliding Drawers for Storage.
- Head Light Band Hanger
- Otoscope Holder
- x2 Endoscope Holders
- Waste Container/Cotton Gauze Box with Lid

Specifications

Light Source:	Illumination: Halogen Power: 15v/150W Color Temperature: 3200K
Suction:	Max Vacuum: 710 ± 10mm Hg regulable through Vacuum Control Knob. Suction Capacity: 40 ltrs. / MIN Vacuum Gauge: Bourden Type 3" 0-760mm Hg calibration Motor: 1 / 4 hp motor
Anti-fog Device:	Airflow Rate: 40 cu ft./min Air Temperature: 55°C
Instrument Trays:	Size: 11x7"
Power Consumption:	500-1000W approx.
Power Supply:	AC 220V-240V, 50/60Hz
Dimensions:	600x580x1335mm (LxWxH) approx.
Weight:	85kg. approx.

ENT OPD UNIT

Model: Basic



Color Options



& more

* Please contact us for more information.

* With a policy of continuous development, we reserve the right to change design and specifications without prior notice.

