

Features

- Fully Flexible Chair with Up-Down, 180 Degree Rotating and Reclining Chair
- German DC CE-Marked Noiseless, Motor-Driven
- Ideally Located Control Panel on the Side of Back Rest
- An Auto-Return Button Returns Patient to the Low, Upright Position
- Foot switch for Up Down Movement
- Long Lasting Hand-sewn Upholstery with High-Density Foam
- Multiple Color Options

Specifications

Height Min / Max	560 / 760mm
Backrest Recline	90 Vertical to 180 deg. Horizontal
Maximum Load	200 kg.
Floor Space	34 x 24"
Rotation	180 deg.
Reclined Position Length	1820mm
Compliance	CE/ ISO 13485
Total Weight	108 kg.
Power Supply	AC 220V-240V, 50/60Hz

ENT Patient Chair

Model: GMENT2000



Color Options   & more

* Please contact us for more information.

* With a policy of continuous development, we reserve the right to change design and specifications without prior notice.



Features

- 360-degree swivel
- Ergonomically Designed Backrest for Special Support for Upper Body & Back
- In-built Gas Cylinder for Effective Seat Adjustment
- Five-castor Base for Smooth Movements
- Long Lasting Hand-sewn Upholstery with High-Density Foam
- Multiple Color Options

Specifications

Height Min / Max	560 / 720mm
Rotation	360 deg.
Seat Diameter	380mm
Back Support	Adjustable Height
Compliance	CE/ ISO 13485
Total Weight	12 kg.

Doctor's Stool

Model: GMDS-2000



Color Options   & more

* Please contact us for more information.

* With a policy of continuous development, we reserve the right to change design and specifications without prior notice.

