


ENT TREATMENT UNIT

Advanced

- Compact Ergonomic Design
- Feature Loaded



GEM SURGITECH



Slim LED X-Ray Viewer
with Adjustable Brightness

2 S.S. Instrument Trays
with Cover

2 Endoscope Holder
Otoscope Holder
Headlight Hanger

Halogen + Digital LED Light Source
with Adjustable Brightness

Anti-Fog Device/Mirror Pre-Heater

Advanced Ear Irrigation System
- Electronic Temperature Control
- Built-in Water Tank
- Adjustable Pressure
- Foot Controlled

Castor Wheels
for Easy Mobility

Standard accessories

- Fibre Optic Cable 4.5mm, 2.3m

ergonomic design

Endoview Endoscopy Camera
with Adjustable Monitor

Spray System

1 Straight + 1 Bent for Liquid Medicine

In-Built Suction Machine

- Twin Jars
- Pressure Gauge & Release Valve

High Performance LED Penlight
for Endoscope Illumination

Thermal Cautery

Sliding Drawers

2 Large Size



Standard colors



French Blue
Barn Red



Rich Black
Gold

Technical Details

Suction:	Max Vacuum: 710 ± 10mm Hg Suction Capacity: 40 ltrs. / Min Vacuum Gauge: Bourden Type 3" 0-760mm Hg calibration Motor: 1 / 4 hp motor Vacuum Pump: Rotary type Jars: Twin Polycarbonate Jars 2+2 ltrs.
Light Source	Outlet A: Illumination: LED Lumens: >200,000 lux Color Temp: 5700K Outlet B: Illumination: Halogen Power: 150W Color Temp: 3200K
Anti Fog Device:	Airflow Rate: 40 cu ft./min Air Temperature: 55° C
Endoscopy Camera <i>Endoview</i>	Image Sensor: 1/3" Sony Sensor Resolution: 1200tvl (720p) Signal: PAL Cable Length: 3 metre
Penlight:	Power: 10W Illumination: LED Color Temp: 6000K
Power Consumption:	500-1500W
Power Supply:	AC 220V-240V, 50/60Hz
Dimensions:	820 x 600 x 1335mm (L x W x H)
Weight:	95kg. approx.

With a policy of continuous development, Gem Surgitech reserves the right to change design and specifications without prior notice.



GEM SURGITECH

Contact

Address: 155, Preet Nagar, Ambala Cantt., Haryana (133001) India

Contact Info: Ph No. +91-999625-6235, +91-98174-19409

Email: info@gemsurgitech.com, gemsurgitech@rediffmail.com

Website: www.gemsurgitech.com

(An ISO 9001:2015, ISO 13485:2016 Certified Company)

Visit our website:

www.gemsurgitech.com