

Endovision Unit 17"

Portable Medical Endoscopy Unit




Screen

- 17" LED Screen
- HD OLED Panel

Case

- Fibre Body
- Solid Structure
- Lightweight



Endoscopy Camera

- EndoView Endoscopy Camera
- Sony Image Sensor
- High Resolution Camera
- Anti-Glare Coated Sensor

Light Source

- 120W Digital LED Light Source
- Lightweight Fibre Optic Light Cable
- Adjustable Brightness



Procedure can be viewed or recorded to external Monitor/TV, Laptop or Smartphone



Windows



Android



Apple

Address: 155, Preet Nagar, Ambala Cantt., Haryana (133001) India

Contact Info: Ph No. +91-9996256235, +91-9817419409

Email: info@gemsurgitech.com, gemsurgitech@rediffmail.com

Website: www.gemsurgitech.com

(An ISO 9001:2015, ISO 13485:2016 Certified Company)


Endovision Unit 17"

Portable Medical Endoscopy Unit



Technical Specifications:

Sensor:	Sony 1/3" Image Sensor
Camera Resolution:	1280 x 720p (HD)
Frame Rate:	60 fps
White Balance:	Yes
Camera Cable Length:	2.3 metre
LED Power:	120W
Brightness:	Adjustable
LED Lux:	>100000lux
Color Temp:	6000K (White)
Screen Size:	17"
Compliance:	CE / ISO 13485
Package Contents:	Main Unit, Camera Head, Fibre Optic Light Cable, Remote, AV Cable, Power Cord, Carry Bag
Power Consumption:	130W
Power Supply:	AC 220V-240V, 50/60Hz
Dimensions:	39 x 33 x 43 cm (W x D x H)
Weight:	6kg.



* Please contact us for more information.

* With a policy of continuous development, we reserve the right to change design and specifications without prior notice.



Address: 155 Preet Nagar, Ambala Cantt., Haryana - 133001
Contact Info: Ph No. +91-999625-6235, +91-98174-19409
Email: info@gemsurgitech.com, gemsurgitech@rediffmail.com
Website: www.gemsurgitech.com