

Endovision Unit

Portable Medical Endoscopy Unit




Screen

13.5" LED Panel
HD Screen

Case

Custom Designed Case
Scratch Resistant and Solid Structure
Lightweight & Lockable



Endoscopy Camera

EndoView Endoscopy Camera
Sony Image Sensor
High Resolution Camera
Anti-Glare Coated Sensor

Light Source

10W LED Light Source
High Performance Bright LED
Adjustable Brightness



Procedure can be viewed or recorded to external
Monitor/TV, Laptop or Smartphone



Windows



Android



Apple

Address: 155, Preet Nagar, Ambala Cantt., Haryana (133001) India

Contact Info: Ph No. +91-9996256235, +91-9817419409

Email: info@gemsurgitech.com, gemsurgitech@rediffmail.com

Website: www.gemsurgitech.com

(An ISO 9001:2015, ISO 13485:2016 Certified Company)

Endovision Unit

Portable Medical Endoscopy Unit



Technical Specifications:

Sensor:	Sony 1/3" Image Sensor
Camera Resolution:	1280 x 720p (HD)
Frame Rate:	60 fps
White Balance:	Preset
Camera Cable Length:	2.3 metre
LED Power:	10W
Brightness:	Adjustable
LED Lux:	>100000lux
Color Temp:	6000K (White)
Screen Size:	13.5"
Compliance:	CE / ISO 13485
Package Contents:	Main Unit, Camera Head, Light Source, Light Cable, Remote, AV Cable, Capture Card, Power Cord, Carry Bag
Power Consumption:	35W
Power Supply:	AC 220V-240V, 50/60Hz
Dimensions:	33 x 30 x 16 cm (L x D x H)
Weight:	5kg.

* Please contact us for more information.

* With a policy of continuous development, we reserve the right to change design and specifications without prior notice.



Address: 155, Preet Nagar, Ambala Cantt., Haryana (133001) India

Contact Info: Ph No. +91-9996256235, +91-9817419409

Email: info@gemsurgitech.com, gemsurgitech@rediffmail.com

Website: www.gemsurgitech.com

(An ISO 9001:2015, ISO 13485:2016 Certified Company)