



Specifications

Model	GMH5	GMH7
LED Power	10W	10W
Battery Capacity	2600mah	2600mah
Working Duration	4 hrs.	4 hrs.
Adjustable Brightness	✓	✓
Color Temp	6000K	6000K
Lux	>100,000Lux	>100,000Lux
Adjustable Light Spot	✗	✓
Extra Battery	Optional	Optional
Magnifying Loupe	Optional	Optional
Package Contents	Headlight, Charger, Dust Cover	Headlight, Charger, Dust Cover
Compliance	CE / ISO 13485	CE / ISO 13485
Charging Input	5V	5V
Weight	190g.	190g.

Features

Lightweight

Soft Cushioning

Bright & Energy Efficient

Quick Charging

* Please contact us for more information.

* With a policy of continuous development, we reserve the right to change design and specifications without prior notice.

