



Specifications

LED Power	150W
Lux	>200,000Lux
Color Temp	6000K
CRI	>95 High CRI
LED Life	60,000 hrs.
Adjustable Brightness	1-100%
Display Size	3.25"
Package Contents	Light Source, Main Power Cord
Fibre Optic Cable	Optional
Compliance	CE / ISO 13485
Power Supply	100-240V AC, 50/60Hz
Dimensions (mm)	220 L x 260 W x 100 H
Weight	2.5 kg.

Features

High Intensity Cold LED Illumination

True Color Rendering (CRI Optimized)

High Luminous Efficiency

Multi-Info Smart Display

* Please contact us for more information.

* With a policy of continuous development, we reserve the right to change design and specifications without prior notice.

