


MEDICAL ENDOSCOPY RECORDER

CAPTURE. RECORD. REVIEW.
WITH CLARITY & RELIABILITY.

- HD**
1080P FULL HD
High quality video recording
- AV & HDMI INPUT**
Supports multiple video sources
- BUILT-IN RECHARGEABLE BATTERY**
Portable & reliable for extended use
- RECORD TO PENDRIVE**
Easy storage & quick data transfer
- 3" DISPLAY**
Clear, compact and easy to use



This Medical Endoscopy Recorder is a compact, *all-in-one* solution designed for medical professionals. It ensures reliable high-definition recording and easy data management, making it ideal for endoscopy procedures.

KEY FEATURES

HD					
1080P FULL HD RECORDING	AV & HDMI INPUT	BUILT-IN BATTERY	RECORD TO PENDRIVE	3" DISPLAY	EASY TO USE
Captures crystal clear videos in 1080P Full HD resolution at 60FPS.	Compatible with a wide range of endoscopes and video sources.	High-capacity rechargeable battery for cordless operation on the go.	Save recordings directly to USB pendrive – no PC required.	3-inch high-resolution display for real-time viewing and easy navigation.	Intuitive controls and user-friendly interface for seamless operation.



Address: 155 Preet Nagar, Ambala Cantt., Haryana - 133001
Contact Info: Ph No. +91-999625-6235, +91-98174-19409
Email: info@gemsurgitech.com, gemsurgitech@rediffmail.com
Website: www.gemsurgitech.com