



Specifications

Max Speed	6000 rpm
Speed	Adjustable
Motor Type	Brushless
Modes	Forward, Reverse, Oscillation
Compatibility	Medtronic & Other Blades / Burrs Compatible
Display Type	Digital
Display Size	3.25"
Package Contents	Handpiece, Main Control Unit, Foot Pedal, Main Power Cord, Blade
Compliance	CE / ISO 13485
Power Supply	100-240V AC, 50/60Hz
Dimensions (mm)	220 L x 260 W x 100 H

Features

Powerful Brushless
Motor

Hi-Speed Oscillation

Digital Display

Multi-Brand Blades
Compatibility

Ergonomic and
Lightweight Handpiece

* Please contact us for more information.

* With a policy of continuous development, we reserve the right to change design and specifications without prior notice.

