



Specifications

Micromotor

Power	45W
Max Speed	35,000 rpm
Torque	2.9 ncm
Package Contents	Control Unit, Handmotor, Foot Pedal, Holder
Compliance	CE / ISO 13485
Power Supply	100-240V AC, 50/60Hz
Dimensions (mm)	118W x 156D x 78H
Weight	1.25kg

Handpiece

Gear Ratio	1:1
Angle	Straight
Max Speed	40,000 rpm
Package Contents	Handpiece
Compliance	CE / ISO 13485
Dimensions (mm)	L 85 x Ø 20
Weight	71 g.

Features

Compact Size

Step-less Speed Control

Forward & Reverse Direction Switch

Low Noise & Vibration

* Please contact us for more information.

* With a policy of continuous development, we reserve the right to change design and specifications without prior notice.

