


# Portable LED Light Source

Model: PLS10+



## Specifications

|                              |                                            |
|------------------------------|--------------------------------------------|
| <b>LED Power</b>             | 10W                                        |
| <b>Adjustable Brightness</b> | 3 Brightness Modes: Low, Medium, Maximum   |
| <b>Lux</b>                   | >150,000Lux                                |
| <b>Color Temp</b>            | 6000K                                      |
| <b>Working Duration</b>      | 4 hrs.                                     |
| <b>Rechargeable</b>          | Yes                                        |
| <b>Extra Battery</b>         | <b>Optional</b>                            |
| <b>Package Contents</b>      | Light Source, Charger, Battery, Carry Case |
| <b>Compliance</b>            | CE/ ISO 13485                              |
| <b>Dimensions</b>            | 25 x 135mm                                 |
| <b>Weight</b>                | 120g. (with battery)                       |

## Features

Ultra Bright LED

Compact, Wireless & Lightweight

Zero Heat Transfer

Super Quick Charging

\* Please contact us for more information.

\* With a policy of continuous development, we reserve the right to change design and specifications without prior notice.

